

Connecting

from

heart

to

heart



SonoTrax Series Ultrasonic Pocket Doppler can meet the routine examination requirements of obstetricians, midwives and pregnant women.

SonoTrax

Ultrasonic Pocket Doppler

Features

- Easy-to-use
- Ergonomic design
- Interchangeable probe
- Built-in audio recorder
- Large backlight LCD display
- Accurate FHR detection with clear sound
- Support 2MHz, 3MHz, 4MHz, 8MHz waterproof probes
- Automatic and manual FHR counting
- Long time continuously use (4 hours)
- Automatic power-off



SonoTrax

Ultrasonic Pocket Doppler



Specifications

Display: 45x25mm LCD display

Heart rate range: 50~210bpm

Resolution: 1bpm

Audio sampling frequency: 4KHz /8KHz

Recording length: 480s/240s

Audio output power: 0.5W

Auto shut-off: 1 minute

Size: 32 (Depth) x 85 (Width) x 138 (Height) mm

Weight: 290±5g (including one battery)

Model	Interchangeable probes	FHR display	Built-in audio recorder	Backlight	Waterproof probe	Probe auto-detection display	Battery status display	Low battery indication	AC recharge
SonoTrax Lite	Yes	No	No	No	Yes	No	No	Yes	No
SonoTrax Basic	Yes	Yes	No	No	Yes	Yes	No	Yes	No
SonoTrax Basic A	Yes	Yes	No	Yes	Yes	Yes	Yes	Yes	No
SonoTrax Pro	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No
SonoTrax II	Yes	Yes	No	Yes	Yes	Yes	Yes	Yes	Yes
SonoTrax II Pro	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Optional accessories	Rechargeable Li-ion battery, charger, earphone, carry bag, 9V alkaline battery, adapter								